EN User Guide

Thank you for purchasing this Lifetrons® product. This user guide outlines the operating process and safety precautions. Please read it carefully to make sure you are using your product properly and retain it for future reference.

How to Use

- A. Ultraschall & Positive Ion Mode**
 Ultrasonic power combined with positive ions uses the characteristics of opposite attraction via ionic charge to export dirt and debris in the pores.
- Use on clean skin and be sure to remove makeup before use.
 - Remove the protective cap ➊. Press and hold the power button ➋ for 2 seconds to turn on the device. The default program will begin with ultrasonic level „L“ for low power and positive ion mode (+) for deep cleansing.
 - Press ➌ to increase the ultrasonic level (see level indicator ➍). Levels : L: Low M: Medium H: High
 - Choose a suitable power level for your skin and start the cleansing process. Spray your face with water. Hold the device with both sides of your hand touching the conductive pole side strips ➎ to activate the positive ion function. Gently move the metal tip of the device ➏ (face down) over your skin to cleanse your face and pores from dirt and debris. Ultrasonic technology combined with positive ions attracts dirt and debris, allowing it to escape from your pores — sometimes causing water atomization. For best results, hold the device face down while sliding it firmly against your skin.

- B. Ultraschall & Negative Ion Mode**
 Ultrasonic power combined with negative ions uses the characteristics of the same ionic charges to impart the nutrition of your skin care products deep into the skin's layers.
- Use on clean skin and be sure to remove makeup before use. Serum, moisturizer or cream (not included) must be applied in combination with the steps below to ensure proper operation and maximum benefits.
 - Remove the protective cap ➊. Press and hold the power button ➋ for 2 seconds to turn on the device. The default program will begin. Press ➌ and hold for 2 seconds to switch to ultrasonic level "L" and negative ion mode ➍.
 - Press ➌ to increase the ultrasonic level (see level indicator ➍). Levels : L: Low M: Medium H: High
 - Apply face cream or serum onto clean skin (not included). Choose a suitable power level for your skin and start the nutrition importation procedure. Hold the tip of the device face up ➎ and hold it with both sides of your hand touching the conductive pole side strips to activate the negative ion function. Gently move the metal tip across your skin for 3 minutes to help impart the nutrition from your skin cream into your pores. For best results, hold the device face up while sliding it firmly against your skin.

C. Electrical Muscle Stimulation (EMS) Mode
 EMS treatment may help tone facial muscles, enhance contours and help skin care nutrition permeate into the underlying skin.

1. Use on clean skin and be sure to remove makeup before use. Serum, moisturizer or cream (not included) must be applied in combination with the steps below to ensure proper operation and maximum benefits.

2. Remove the protective cap ➊. Press and hold the power button ➋ for 2 seconds to turn on the device. The default program will begin. Press and hold "EMS" for 2 seconds to switch to the low level "L" EMS mode.

Press "EMS" to increase the EMS level (see level indicator ➍). Levels : L: Low M: Medium H: High

4. Choose a suitable power level for your skin and start the EMS treatment. Turn the tip of the device face up to begin the EMS treatment. Starting at the bottom of the face, move the tip in a continuous upward direction on the skin ➎.

D. How to Shut Down
 1. The product will automatically shut down after 15 minutes of continuous use. It will beep 14 times (once every minute) to help you keep track of the usage. If you want to shorten the usage time, simply press and hold the power button for 3 seconds. The device will beep 3 times before shut down.

E. How to Recharge
 1. Open the USB cover.
 2. Connect ➉ (Micro-USB jack) to "D3" (Micro-USB port).
 3. Connect ➊ (USB jack) to a USB port with a power supply (e.g., AC adaptor with USB power supply).
 4. A quick flashing blue LED light will appear when the device is charging. Once the charge cycle is complete, the solid blue light will disappear and you can disconnect the device.

F. LED Indicator
Low battery: Flashing blue LED light every 5 seconds
Charging: Flashing blue LED light appears every second
Fully charged battery: Solid blue LED light appears

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