

How to Use

A. Ultrasonic Mode

Ultrasonic power uses the characteristics of opposite attraction via the ionic charge to export dirt and debris in the pores.

1. Use on clean skin and be sure to remove makeup before use.
2. Press the power button ❶ to turn on the device. Press ❸ to start the ultrasonic mode.
3. Gently press the massage head ❷ onto your skin to start the facial massage. Slowly and gently glide the device over your entire face with light pressure.
4. Keep pressing ❹ to choose your desired ultrasonic level (see level indicators ❹).
Levels
L: Low
M: Medium
H: High

Note: Remove any makeup before use. Choose a suitable power level for your skin and start the cleansing process. Spray some water on your face. Do not treat the eye area with ultrasonic mode.

B. Positive Ion Mode

1. Use on clean skin and be sure to remove makeup before use.
2. Press the power button ❶ to turn on the device. Press ❺ to start the positive ion mode, red LEDs appear ❹.
3. Hold the device with both sides of your hand touching the conductive pole side strips ❻ to activate positive ion mode.

C. Negative Ion Mode

1. Use on clean skin and be sure to remove makeup before use.
2. Press the power button ❶ to turn on the device. Press ❻ 2 times to start the positive ion mode, green LEDs appear ❹.
3. Hold the device with both sides of your hand touching the conductive pole side strips ❻ to activate negative ion mode.

D. How to Shut Down

1. The product will automatically shut down after 15 minutes of continuous use. It will beep 14 times (once every minute) to help you keep track of the usage time. If you want to shorten the usage time, simply press and hold the power button. The device will beep 3 times before shutting down.

E. How to Recharge

1. Plug the battery charger into the charging port ❸.
2. Connect the charger (e.g., AC adaptor) with a power supply

