CM-300

EN User Guide

Thank you for purchasing this Lifetrons® product. This user guide outlines the operating procedures and safety precautions. Please read it carefully to ensure you are using your product properly and retain it for future reference.

How to Use

Always use this product on clean, dry skin. Start with the exfoliating treatment, then use the blue light and finally the facial massage.

A. Changing the Treatment Heads

The treatment heads attach magnetically to the top of the device. Do not change the treatment head when the power is on. Remove any makeup and dry your skin before use. Do not use on wet skin.

1. There are 3 treatment heads included with this device: normal exfoliating head, sensitive exfoliating head and blue light head.

2. Select the treatment head you would like to use and attach it to the body of the device.

B. Using the Exfoliating Heads

1. Attach the exfoliating head of your choice to the body of the device. Press the power button. 1 time to turn on the device. The device will automatically shut down after 10 minutes of continuous use.

2. Gently move the exfoliating head on your face until you feel a slight tingle on your skin. Never use too much force and always avoid the eye area. Make sure that there is always full skin contact between the exfoliating head and the skin.

3. Use the hand to hold the skin taut while moving the device across your skin with the other hand. This makes it easier to achieve proper suction, especially on the nose and cheek area. Gently move the exfoliating head in straight lines upwards. Start at the chin area, then move to the cheeks and forehead. Note: Depending on your skin type, you may experience a warm and tingling sensation with mild skin redness. These effects are temporary and usually diminish within 5-10 minutes.

4. Treatment heads can be cleaned by wiping them. We recommend treating each head 2-3 times before the first treatment. After the second treatment, you can increase treating each area 4-6 times or start using the normal exfoliating head. We recommend using the exfoliating function bi-weekly. More frequent treatment may cause skin irritation.

Note: We recommend using the sensitive treatment head first to see how your skin reacts. You can do a test on the back of your hand before applying it to your face. If your skin reacts well to the sensitive head, you can try using the normal treatment head in the future.

C. After the Exfoliation Treatment

Apply a moisturizing mask or face cream after the treatment and always use sun protection when outdoors. Do not use any other method of mechanical exfoliation, chemical exfoliation, under-eye electrolysis, facial waxing or use depilatories within 3 days after the treatment.

D. Using the Blue Photon Head

1. Place the blue photon light head on the body of the device. Press the power button. 1 time to turn on the device. Then use the suction when the blue light is on. The device will automatically shut down after 10 minutes of continuous use.

2. Place the blue photon light head in direct contact with the skin wherever you want to treat acne, blackheads, oil or blemishes on your skin.

3. The device will automatically stop massaging when it is no longer in contact with the skin. After use, gently clean the massage head with a cotton pad.

Note: Skin must be clean before use. Serum, moisturizer or cream (not included) must be applied in combination with the negative ion facial massager to ensure proper operation and maximum benefits.

E. Cleaning the Treatment Heads

1. Remove the treatment head before cleaning. Rinse the device with warm water to remove any residue and let them thoroughly dry before using again. Wipe the blue photon light head with a clean dry cloth. Do not use any liquid to clean the device body. Keep in mind that the device should only be wiped with a clean dry cloth. Never use cleaning pads, abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean this device.

F. Using the Facial Massager

1. Always use clean skin. Apply a thin layer of serum, moisturizer or cream to your face and eye area.

2. Hold the device with both sides of your hand touching the conductive pulse strips (1) and the conductive needels (2) on your skin to start the facial massage. Slowly and gently move the device on your face covering all areas with light pressure.

3. The device will automatically stop massaging when it is no longer in contact with the skin. After use, gently clean the massage head with a cotton pad.

Note: Serum must be clean before use. Serum, moisturizer or cream (not included) must be applied in combination with the negative ion facial massager to ensure proper operation and maximum benefits.

G. Recharging the Device

1. Open the micro-USB charging port cover. Connect part ① (micro-USB) to part ② (micro-USB port).

2. Connect part ③ (USB jack) to a USB port with power supply (e.g., AC adapter with USB power supply).

H. Battery Indicator

Low battery: Red LED lights three times and the device turns off (only visible when power is on)

When charging: Green LED lights once per second

Fully charged battery: Solid green LED lights appear

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